

Inogen One G2 Batterie Laufzeiten

12 Cell-Batterie

Setting 1.0	4.0 Stunden
Setting 2.0	3.5 Stunden
Setting 3.0	2,7 Stunden
Setting 4.0	2.2 Stunden
Setting 5.0	2.0 Stunden

24 Cell-Batterie

Setting 1.0	8.0 Stunden
Setting 2.0	7.0 Stunden
Setting 3.0	5.4 Stunden
Setting 4.0	4.4 Stunden
Setting 5.0	4.0 Stunden

Inogen One G3 Batterie Laufzeiten

8 Cell-Batterie

Setting 1.0	4.5 Stunden
Setting 2.0	4.0 Stunden
Setting 3.0	3.0 Stunden
Setting 4.0	2.0 Stunden

16 Cell-Batterie

Setting 1.0	9.0 Stunden
Setting 2.0	8.0 Stunden
Setting 3.0	6.0 Stunden
Setting 4.0	4.0 Stunden